
CAROUSEL Temperature Controller

RD RADLEYS
DISCOVERY
TECHNOLOGIES

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 61 010-1; EN 50 081; EN 50 082; EN 55 011 and EN 60 555.

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
Garantee

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer.

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Carousel Temperature Controllers are manufactured exclusive for Radleys Discovery Technologies Ltd by IKA LABORTECHNIK Janke & Kunkel GmbH & CO. KG

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Contents

	Page
Correct use	8
Unpacking	8
Commissioning	8
System configuration	8
Operating modes	9
Setting the nominal temperature	9
Setting the HI-TEMP	9
Replacing extension cable and sensor	10
Error messages	10
Analog output	10
Accessories	10
Technical data	10
Spare parts list	14
Spare parts diagram	15

Correct use

The CAROUSEL Temperature Controller is an instrument for the precise control of temperatures. The sensor of the instrument must be immersed in at least 20mm deep in the medium to be tempered.

The instrument can be connected to any magnetic agitator with contact thermometer connection which conforms with the requirements listed in the Technical Data (see section Technical Data).

Unpacking

Please unpack the machine carefully and inspect for damage. It is important that any transit damage should be noted at the time of unpacking. In certain circumstances it may be necessary to investigate immediately (post, rail or freight forwarder).

As supplied, the package contains: one CAROUSEL Temperature Controller and one set of operating instructions.

Commissioning

Plug the instrument into the contact thermometer jack on the agitator (e.g. CAROUSEL Stirring Hotplate). The instrument is activated when the agitator is switched on. If there is no voltage supply between 8V and 16V DC at the contact thermometer jack, the instrument must be run off the H52 power supply kit (H50 twin-jack adapter and H51 plug power supply). The instrument will then receive the required voltage from the power supply or twin-jack adapter.

An automatic self-test is performed when the instrument is activated. During the self-test all the segments of the LCD display are visible (**Fig. 1**). Then the LCD displays illustrated in **Fig. 2** to **Fig. 7** appear (standard factory setting).

System configuration

Before starting to work with the CAROUSEL Temperature Controller the system configuration must be set. To do this, hold down the "-" key when switching on.

1. Setting the type of sensor used

The LCD display indicates the set type of sensor, e.g. H 62 (**Fig. 8**). Accessories see page 10. Use the "+" and "-" keys to set the type of sensor H 60 ... H 67 you are using. To change the setting, hold down one of the two keys for approx. 2 seconds. The displayed sensor is set and stored if no key is pressed for approx. 5 seconds.

2. Setting the extension cable used

The LCD display indicates the set type of extension cable (**Fig. 9**) in small letters. The setting is made by pressing the "+" and "-" keys. Possible settings are: H 70, H 71 or — if no extension cable is used. The displayed extension cable is

set and stored if no key is pressed for approx. 5 seconds.

3. Setting the control mode

The LCD display indicates the control setting. Three control modes are available. The setting can be changed by pressing the "+" and "-" keys. The displayed control mode is set and stored if no key is pressed for approx. 5 seconds.

- a) Fuzzy control (**Fig.10**). Standard factory setting, best control behavior.
- b) Fast Fuzzy control (**Fig.11**). Faster warming up than with a), very good control behavior in stabilized state. ">" displays Fast Fuzzy control mode.
- c) Two-Point control (**Fig.12**). Faster warming up than with a) and b), severe overshwing on first warm-up.

Possible application: Using the CAROUSEL Temperature Controller with KSR 11 as a limit monitor. The heater or instrument remains switched on until the set nominal temperature is exceeded. "" signifies Two-Point control mode.*

Operating modes

The CAROUSEL Temperature Controller can be run in three different operating modes. (Displays see **Fig.13** to **Fig.15**)

Operating mode "A" The instrument is factory-set to this mode. When the instrument is switched on the nominal temperature is always set to -10°C. The nominal temperature can be adjusted with the "+" and "-" keys. The HI-TEMP (maximum nominal temperature setting) can only be set in this mode.

Operating mode "B" In this mode the most recently set nominal temperature is adopted and becomes active when the instrument is switched on again. The setting can be adjusted with the "+" and "-" keys.

Operating mode "C" In this mode the current settings are adopted and become active when the instrument is switched on again. The setting cannot be adjusted. Holding down the "+" key when switching on the instrument switches to the next operating mode in the sequence A-B-C-A-B-C-A etc.

Setting the HI-TEMP

The "HI-TEMP" (maximum nominal temperature setting) can only be set in operating mode "A". After the self-test, HI-TEMP appears in the LCD display (**Fig.16**). The safe temperature can now be set between 0 and 400°C using the "+" and "-" keys. The displayed HI-TEMP is set and stored if no key is pressed for approx. 5 seconds. The LCD display then shows "HI OK" (**Fig.17**).

Setting the nominal temperature

The nominal temperature is adjusted in 1K steps by pressing the "+" and "-" keys. If the key is held down for longer than 5 seconds, the nominal temperature setting changes in 10K steps. The LCD display shows: the current actual temperature with physical units, the nominal temperature at top right, and OFF to indicate that the heater of the connected device is switched off. (**Fig.18**)

If the set nominal temperature is higher than the current actual temperature, the word OFF in the display changes to ON and the heater of the connected device is switched on. (**Fig.19**)

The star (*) in the LCD display indicates that the instrument is being operated in mode "B" or "C". (**Fig.20**)

Replacing extension cable and sensor

In order to replace the sensor or the extension cable, pull the sensor and the cap downward until the plastic snap fastener is released. Then push the sensor or extension cable with the cap over the fitting on the instrument to reestablish the connection.

Error messages

The following error messages may occur (**Fig.21**) and have the causes listed below:

- Er 1 Sensor not connected or defective (open circuit).
- Er 2 Analog output loaded with <math><50\text{k}\Omega</math>, output voltage too low.
- Er 3 Internal temperature above $+76\text{ }^\circ\text{C}$.
- Er 5 Sensor not in medium (no change in medium temperature after heater has been on for 3 minutes).
- Er 6 Sensor short circuit.
- Er 7 The maximum set SAFE TEMP is exceeded.
- Er 8 Calibration error (re-calibration by field service required).
- Er 9 IC107 (4066) defective
- Er 10 T₁₀₅ defective, OUT_TEMP $>480\text{ }^\circ\text{C}$ (4.8V)
- Er 11 Supply voltage $<7.5\text{V}$
- Er 12 Supply voltage $>16.1\text{V}$

Analog output

The actual temperature is represented as a voltage (1V corresponds to 100°C) at the 3.5mm jack of the twin-jack adapter.

Tolerance: $0\text{...}100\text{ }^\circ\text{C} \pm 10\text{ mV}$
 $100\text{...}400\text{ }^\circ\text{C} \pm 25\text{ mV}$

The internal resistance of the connected measuring instrument or recorder must be greater than $50\text{ k}\Omega$.

Accessories

RR98107	Holding Rod
RR98106	Boss Head
-	H 50 Twin-jack adapter
-	H 52 Power supply kit
RR98100	H 60 Short sensor 60 mm stainless steel
RR98098	H 62 Standard sensor 270 mm stainless steel
-	H 64 Glass sensor 50 mm glass
RR98099	H 66 Glass sensor 270 mm glass
RR98102	H 70 Extension cable, 1 m long
-	H 71 Extension cable, 2.5 m long
-	Ak 2.5 Analogue cable

Technical data

Measuring range:	$^\circ\text{C}$	-10 ... 400
Resolution:	K	0.1 (-10 ... +100 $^\circ\text{C}$)
	K	1 (100 ... 400 $^\circ\text{C}$)
Measuring accuracy:	K	± 0.3 + sensor tolerance Pt 1000 DIN IEC 751 Class B
Setting accuracy:	K	1
Control deviation:	K	± 1
Supply voltage:	V DC	8 ... 16
Current consumption:	mA	16 (at 9V)
Ambient temperature:	$^\circ\text{C}$	0 ... 60
Relative humidity:	%	80
ON time:	%	100
Degree of protection:		IP 44
Cable length:	mm	350
Plug:		6-pin DIN 45 322
Connection:		DIN 12 878 Class 2
Dimensions: (W x D x H)	mm	70 x 150 x 65 (without sensor)
Weight:	kg	0.2

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

List of spare parts

Item	Designation
1	Rear wall
2	Clamp
3	Sealing ring
4	Spiral cable
5	Sealing insert
6	Union nut
7	Housing
8	Dished washer
9	Membrane keypad
10	Sensor jack; complete
12	BLP temperature regulator
13	Plastite screw 2.5 x 8
14	H 62 stainless steel sensor Pt1000
15	Cap
16	Gasket

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
Carousel Stirring Hotplates and Accessories

RR98072	Carousel Stirring Hotplate 230 volt & UK Plug
RR98072/EURO	Carousel Stirring Hotplate 230 volt & Euro Plug
RR98072/SWISS	Carousel Stirring Hotplate 240 volt & Swiss Plug
RR98072/USA	Carousel Stirring Hotplate 115 volt & US Plug
RR98072/JAP	Carousel Stirring Hotplate 110 volt & Japanese Plug
RR98073	Fuzzy Logic Digital Temperature Controller
RR98098	Spare Stainless Steel Sensor for Temperature Controller
RR98099	Glass Coated Sensor for Fuzzy Logic Temperature Controller
RR98075	Small Cross Shaped Stirring Bar, Pk 40
RR98091	Rare Earth Medium Cross Shaped Stirring Bar, Pk 20
RR98070	Large Octagonal Stirring Bar, Pk 20
RR98071	Small Octagonal Stirring Bar, Pk 40
RR98096	Rare Earth Small Elliptical Stirring Bar, Pk 40
RR98097	Rare Earth Large Elliptical Stirring Bar, Pk 20
RR98095	PTFE Stirring Bar Evaluation Kit, Box 18
RR98103	Stainless Steel Bath with Handle
RR98104	Stainless Steel Bath
RR98105	Special Sand Bath Filling
RR98094	Stirring Bar Retriever, 350mm
RR98106	Boss Head
RR98107	Holding Rod
RR98108	Support Rod
RR98109	Side Support Bracket
RR98110	Holding Fork and clamp
RR98111	Supporting Plate, 200mm diameter
RR98112	Thermomag Circulator
RR98100	Short Stainless Steel Sensor for Fuzzy Logic Temperature Controller
RR98102	1 metre Extension Cable for remote temperature sensing

**RADLEYS
DISCOVERY
TECHNOLOGIES**

Radleys Discovery Technologies Ltd
Shire Hill, Saffron Walden,
Essex, CB11 3AZ. United Kingdom.
Tel: +44-(0)1799-513320
Fax: +44-(0)1799-513283
Email: sales@radleys.co.uk
Web: www.radleys.com