Elektronisches Kontaktthermometer
Electronic Contact Thermometer

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
Please read this manual carefully before using your new electronic contact thermometer.

The operating manual provides clear and simple instructions for using the thermometer.

Sections that are useful and important for understanding how the device works are highlighted by bars in the instruction text.

For safe handling of the thermometer, please follow the safety notes highlighted with the symbol .

Safety Notes

Proceed with caution when measuring high temperatures! Do not touch the temperature sensor with bare hands immediately after measuring! Danger of burn injuries! Never put the thermometer in the hot zone during measurement. Risk of explosion of the lithium battery!

Do not use the thermometer in explosion hazard areas! Risk of accident in case of an explosion!

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
Table of Contents

1	Description	7
2	Scope of Delivery	9
3	Installation	11
4	Operation	13
4.1	Setting the target temperature ..	13
4.2	Display of maximum temperature	15
4.3	Switching off the contact thermometer	15
5	Changing the Sensor	17
6	Troubleshooting	19
7	Calibration Service.....	21
Appendix		21
I	Accessories	21
II	Technical Data	23

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1 Description

The electronic contact thermometer GFX 393 / GFX 393 B is used for temperature adjustment with heating plates and similar laboratory equipment.

Its microprocessor-controlled electronics ensure high measuring precision and linearity over the entire measuring range. Fuzzy logic control guarantees optimal heating with low overshoot.

The unit features a safety and control circuit in accordance with DIN 12878 classes 1 and 2, as well as sensor short circuit and breakage protection.

- 1 Display (LCD)
- 2 Control-LED
- 3 MAX key
- 4 SET key
- 5 Temperature sensor
- 6 Key  (higher target temperature)
- 7 Key  (lower target temperature)
- 8 5-pole DIN plug

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- Netzadapter H 37
- Relaisbox RB 1691
- Zubehör

2 Scope of Delivery

Check whether the contents of the packaging are complete and undamaged. If you discover any damage or have any reason for complaint, please contact:

ebro Electronic GmbH & Co. KG

Peringerstr. 10

85055 Ingolstadt

Phone: +49-(0)841-9 54 78-0

Fax: +49-(0)841- 9 54 78 80

E-mail: info@ebro.de

Internet: <http://www.ebro.de>

Parts List:

- Contact thermometer GFX 393 with calibration certificate

or

- Contact thermometer GFX 393 B with calibration certificate

Depending on order and application:

- Mains adaptor H 37
- Relay box RB RB 1691
- Accessories

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
3 Installation

The power supply required for your GFX 393/393 B is 9 ... 12 V stabilised DC.

The connections of the 5-pole DIN plug are arranged as follows (see illustration):

Insert the DIN plug into the appropriate socket of the laboratory equipment.

The GFX 393/393B is designed for direct connection to IKA heating plates.

For use with heating plates manufactured by Heidolph, the mains adaptor H 37 is available as an optional accessory. For connecting the thermometer to heating plates or other heaters without a DIN socket, the relay box RB 1691 will be required.

The contact thermometer is immediately ready for use and will perform a system test. During this test, all segments of the display will be visible.

Subsequently the measuring channel is tested. The display will show "38.8".

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
After this, the first measured value will appear in °C. The control LED will light up red or green.

If an error message is displayed, please refer to chapter 6, "Troubleshooting".

Position the sensor.

To avoid faulty measurements, use the sensor only as directed.

4 Operation

4.1 Setting the target temperature

Press the »SET« key. The control LED will light up green.

The display will show the most recently set target temperature and two triangles.

Now you can use the arrow keys to adjust the target temperature as required.

Press key briefly:

Value will be increased or decreased in increments of 1° C.

Hold down key:

Value will be increased or decreased in increments of 10°C.

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If no key is pressed, the thermometer will automatically return to measuring mode after approx. 15 seconds.

Control LED red --> heating on

Control LED green --> heating off

4.2 Display of maximum temperature

The maximum temperature reached since switching on the contact thermometer will be saved automatically.

Press the »MAX« key. The display will show the maximum temperature and two triangles pointing down.

To return to measuring mode, press the »MAX« key again.

To cancel the saved value,

- press the »SET« key or,
- turn off the power supply.

4.3 Switching off the contact thermometer

For switching off the contact thermometer, the power supply must be turned off.

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5 Changing the sensor

The thermometer has been calibrated by the manufacturer to the sensor that is supplied with it. Using a different sensor or using an extension cable with the thermometer may result in reduced measurement accuracy.

To unlock the connection, use the supplied open-end spanner to turn the safety cap slightly to the left.

Pull down the safety cap. Now you can see the plug-in connector for the sensor.

Pull the sensor out downwards.

Plug in a new sensor or an extension cable, if required.

Ensure a tight connection. The plug must snap in.

Push the safety cap over the plug-in connector and turn it to the right to lock the connection.

If required, connect the extension cable to the temperature sensor in the same manner, using the safety cap provided for this purpose.

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
6 Troubleshooting

Possible cause	Remedy
Measuring range exceeded	Observe measuring range
Plug-in connection socket – plug	Check connection; plug has to snap in
Sensor defective	Return unit to manufacturer for repair
Temperature below measuring range	Observe measuring range
Sensor short-circuit	Return unit to manufacturer for repair

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7 Calibration Service

In order to guarantee a high measurement accuracy, the thermometer must be calibrated once per year. For this purpose, **ebro Electronic GmbH & Co. KG** offers a calibration service.

- Please fill in the enclosed service coupon.
- After one year we will recall your thermometer for calibration.
- The calibrated thermometer will be returned to you within one week..

Appendix

I Accessoires

Article name	Type
Support	FX 10-2
Sensor clamp	HK
Relay box	RB 1691
Mains adaptor	H 37
Sensor, L = 225 mm	•
Stainless steel, Ø 3 mm	GTF 250 E-FX
Glass-coated, Ø 8 mm	GTF 250 G-FX
Extension cable	
PVC, 1 m	EB 50 K-FX
Teflon, 1 m	EB 53 K-FX

Probes and cables in various lengths available on request.

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* When extension cables are used, protection of the equipment against splashed water is not guaranteed

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