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


## VarioPrep columns

### Instructions for the operation of VarioPrep columns for preparative HPLC

#### Description of the VarioPrep system:

VarioPrep columns feature the possibility to axially adjust both end parts. These end parts are designed as follows:


When designing these columns special care was taken to use only highest quality and completely inert materials. For this reason the sealing combination is made from PTFE for the supporting ring and from the perfluorinated MN-inert elastomer for the sealing ring.

Sealing of the column is achieved by prestretching the elastomer ring; consequently the force, with which the nut is tightened, does not have any effect on the tightness of the system. However, this force has a direct influence on the pressure, with which the packing material is compressed in the column. Thus the nuts are only used for retaining the HPLC material between the piston fittings: if it is tightened too strongly, the HPLC packing can be damaged. Due to the design of the end parts the piston fitting has to be held with an open end wrench when connecting the capillary fittings in order to prevent turning.

For shifting the piston fitting, the capillary connections have to be removed. If this is not done, overpressure or underpressure is formed in the column depending on the direction in which the piston is shifted. Rapid equalisation of this pressure difference can also cause damage to the packing.

Procedure for adjusting a piston fitting in order to remove a dead volume in the column:

VP e 03/07 A028418

After the careful packing and individual testing (see test certificate) at MACHEREY-NAGEL, VarioPrep HPLC columns are supplied with zero dead volume at the column inlet (fig. A). During operation of the column dead volume (fig. B) can be formed as a consequence of the operating conditions, independent of the quality of the silica and the quality of the initial packing.

Dead volume can be repaired as follows (fig. C):

1. Test the column with a standard test mixture (e.g. the so-called aromatics test).
2. Turn off the HPLC pump if the column has to be pressure-free.
3. Remove the connecting capillary. Hold the piston fitting with an open end wrench in order to prevent turning.
4. Loosen the locking screw of the nut.
5. Turn the nut clockwise by hand, until you feel a slight resistance. Please, take care not to turn too far!
6. Tighten the locking screw and connect the capillary again. Hold the piston fitting with an open end wrench in order to prevent turning.
7. Check the quality of the column with the above-mentioned standard test mixture:

if you are satisfied → you can proceed with your work

if you are not → repeat the above procedure, until you get required results (take care! do not tighten the nut too strongly) or until you reach the end of the adjustment length of your column.  
MACHEREY-NAGEL offers interesting conditions for the REFILL of these columns!

Some instructions for opening the VarioPrep columns:

1. Test your column with a standard test mixture prior to opening.
2. Remove the capillary connection to avoid formation of overpressure or underpressure during movement of the piston fitting.
3. Move the piston fitting only carefully: if it is tight, fill the space between the piston fitting and the column wall with, e.g. methanol. Pull the piston fitting, while turning it at the same time (this will result in a screw-shaped movement).
4. If you wish to replace the phase sealing elements (718850 for 10 mm ID columns, 718872 for 21 mm ID columns), please take care not to damage the packing.
5. If you have problems with the tightness of the column, you need to replace the MN-inert sealing combination (718848 for 10 mm ID columns or 718870 for 21 mm ID columns), it should be wetted with methanol before. Relative to the PTFE supporting ring, the MN-inert sealing ring has to be located towards the middle of the column. It is important to take care, that the MN-inert sealing ring is not damaged, when the piston fitting is introduced! If after this replacement, the VarioPrep column is still not yet tight, please contact MACHEREY-NAGEL for further help.  
Please take care, that the surface of the column packing is completely level, before you insert the new phase sealing elements. Of course, the surface has to be perpendicular to the column axis.
6. Tighten the nut by hand, until you feel a clear resistance.
7. After these procedures, test the column under the same conditions as above (step 1).

MACHEREY-NAGEL

MACHEREY-NAGEL

Germany and International: MACHEREY-NAGEL GmbH & Co. KG · Neumann-Neander-Str. 6-8  
D-52355 Dueren · Tel. 02421/969-0 · Fax. 02421/969-199 · e-mail sales@mn-net.com

USA: MACHEREY-NAGEL Inc. · 2850 Emrick Boulevard · Bethlehem, PA 18020  
Tel. 484-821-0984 · Fax 484-821-1272 · e-mail sales-us@mn-net.com

Internet: [www.mn-net.com](http://www.mn-net.com)