
CAROUSEL Stirring Hotplate

RD RADLEYS
DISCOVERY
TECHNOLOGIES

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CE-DECLARATION OF CONFORMITY

GB

We declare under our sole responsibility that this product corresponds to the regulations 89/336EEC and 73/023EEC and conforms with the standards or standardized documents EN 60 010-1; EN 50 082; EN 55 014; EN 55 114 and EN 60 555.

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
Garantee

In accordance with IKA guarantee conditions, the guarantee period is 24 months. For claims under the guarantee please contact your local dealer.

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Carousel Stirring Hotplates are manufactured exclusive for Radleys Discovery Technologies Ltd by IKA LABORTECHNIK Janke & Kunkel GmbH & CO. KG

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Contents

	Page
Guarantee	2
Safety instructions	8
Proper use	8
Unpacking	9
Commissioning	9
Starting	9
Heating Function	9
Controlling the Medium Temperature via Carousel Temperature Controller or Contact Thermometer	10
Stirring Function	11
Maintenance and Cleaning	11
Associated standards and regulations	12
Technical data	12
List of spare parts Carousel Stirring Hotplate	18
Spare parts diagram Carousel Stirring Hotplate	19
Accessories	20

Safety instructions

Caution: Never heat up any medium with ignition point higher than the adjusted safety temperature limit. Value for safety temperature limit see chapter "Technical data".

The instrument may not be operated in rooms with explosion hazard.

Furthermore, it has to be placed upon a fire-proof and/or non burnable surface.

Please make sure that the mains cable does not contact the heating plate!

Attention - Magnetism! Effects of the magnetic field have to be taken into account (e.g. data carriers, cardiac pacemakers).

When using PTFE-coated magnetic bars, the following has to be noted:

Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline - earth metals, as well as with fine-particled powders of metals of the 2. and 3. group of the periodical system at temperatures above 300-400°C.

Only elementary fluorine, chlorine trifluoride und alkaline metals do attack PTFE, halogen hydrocarbons have a reversibly swelling effect. Source: Römpps Chemie-Lexikon and "Ullmann" Bd.19

Exercise caution when touching the housing parts and the heating plate!

Risk of burns! The heating plate can reach temperatures in excess of 300°C.

Please consider possible endangering due to release of toxic or flammable gases caused by heated medium.

Proper use

The Carousel Stirring Hotplate is a stirring instrument with heating function. It is used in laboratories of the pharmaceutical and chemical industry, as well as in universities and schools. The instrument is suitable for tempering substances which are placed on the heating plate in vessels. By means of a built-in stirring drive and with the help of a magnetic bar within the vessel, the substances can be stirred at the same time. The mixing intensity is dependent on the motor speed and on the size of the magnetic bar.

Unpacking


Please unpack the equipment carefully and check for any damage. It is important that any damage which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be set immediately (to your supplier).

The delivery scope covers: a Carousel Stirring Hotplate and operating instructions.

Starting

After the instrument has been switched on, it is in "standby"-operation. The functions "heating" and "stirring" are started with the help of switches on the left and right bottom of the housing

(see illustration 1)

Commissioning

Check whether the voltage specified on the type plate matches the mains voltage available. The power socket used must be earthed (protective earth conductor contact). If these conditions are met, the device is ready to operate after plugging in to the mains plug. If these procedures are not followed, safe operation cannot be guaranteed and/or the equipment may be damaged.

Observe the ambient conditions (temperature, humidity, etc.) listed under Technical Data.

Heating Function

The temperature of the heating plate of the instrument is kept constant by a control circuit. In addition, it is monitored by a safety circuit. The required temperature sensor - a thermocouple - is installed in the heating plate.

The instrument comprises of a silumin heating plate with a heating power of 600 watts. The surface of the heating plate is polished in order to facilitate cleaning. This surface gets coated rather quickly with an aluminium oxide protective layer which is chemically very resistant.

The function "heating" is switched on and/or off by its corresponding switch. When the function "heating" is switched on, the LED "temp" shines green. The temperature of the heating plate is set on the rotary knob "temp" with its corresponding dial (0-300°C). The heating plate is heated up to the temperature set on the rotary knob. During the heating phase, the colour of the LED "temp" changes from green to orange (orange = supply of energy to heating plate; green = set temperature is reached.).


In the event of a disturbance in the control circuit, the heating plate is **permanently cut off** by the safety circuit.

PLEASE NOTE!

The heating plate only heats up, if the contact plug has been plugged (the contact plug is undetachably fitted to the instrument), or if a **Carousel Temperature Controller** and/or a **contact thermometer** has been connected.

(see illustration 2)

illustration 2

Controlling the Medium Temperature via the Carousel Temperature Controller or a Contact Thermometer

The control of the medium temperature via the Carousel Temperature Controller or contact-thermometer should be preferred. Its advantages are: just a short heating-up phase after setting the rated temperature, practically no temperature drift and only a minor temperature waviness. The medium temperature is optimized with the aid of a **FUZZY Logic Controller**, that means, that in steady state any overshooting of the temperature is reduced to a minimum.

The rear of the instrument houses the diode bushing for connection to the Carousel Temperature Controller or contact thermometer, or for plugging the contact plug (see illustration 2).

The electronics of the instrument supply a test current that has to flow over plug pins 3 and 5 of the diode bushing in order for the heating plate to heat.

Safety contact thermometers

acc. to DIN 12 878 class 2 or acc. Gerstel are connected with a three-wire cable, the test current flows through the contact thermometer.

Safety function:

If the test current is interrupted because of e.g. breakage of contact thermometer or falling out of the cable plug, the heating cuts off.

The required plug allocation is indicated on a reference plate above the diode bushing.

Contact thermometer without safety circuit

acc. to DIN 12 878 class 0. The instrument only heats if the test current circuit is closed by an electrical connection of the plug pins 3 and 5.

2-wire connecting cables:

Connect plug pins 3 and 5 of the instrument plug.

3-wire connecting cables:

Here the test current circuit can also be produced in the terminal head of the contact thermometer (connect plug pins 2 and 3).

A 3-wire cable with the required bridge is available (accessories).

The **rated temperature** is set on the **contact thermometer** by means of the rotary magnetic head, whereby the upper edge of the travelling nut is the reference point on the rated value scale.

The required **surface temperature of the heating plate** is pre-selected on the **temp-rotary knob** of the instrument.

If the temp-rotary knob is set to the maximally adjustable temperature, the heating-up takes place as rapidly as possible, yet the medium temperature overshoots the rated value set on the contact thermometer.

If the temp-rotary knob is set to twice the rated value (e.g. with a rated value of 60 °C the temp-rotary knob is set to 120 °C), the result is a good compromise between rapid heating-up and overshooting.

If the temp-rotary knob is set precisely to the rated temperature, the medium does not reach the rated temperature owing to a drop in temperature between heating plate and medium.

On the **Carousel Temperature Controller** the **rated temperature** is set by pressing the rotary knob and by turning it at the same time. When using the ETS-D, there is no overshooting of the medium temperature owing to an "aiming-off allowance" technique.

Stirring Function

The stirring function is started and stopped by its corresponding switch. The motor speed is set on the rotary knob "stirring" with its dial from 1 to 10. If the motor is switched on, it runs-up to the speed set on the rotary knob, and the LED "stirrer" beside the imprint "stirring" shines green.

The dial from 1 - 10 conforms with a speed range from 50-1100 1/min. When idling, the speed is 1100 1/min at maximum setting of the rotary knob (right-hand stop).

If the viscosity increases during the process, the speed drops a little. If the viscosity decreases, e.g. due to adding thinly fluid media, the speed increases a little. By a corresponding resetting of the rotary knob, both of these cases can be balanced, as long as the torque does not exceed its maximum.

Please note that mains fluctuations within the permissible tolerance also effect minor fluctuations of the speed.

Maintenance and Cleaning

The CAROUSEL Stirring Hotplate is maintenance-free. It is subject only to the natural wear and tear of components and their statistical failure rate.

When ordering spare parts, please give the manufacturing number shown on the type plate, the machine type and the name of the spare part.

Please send in equipment for repair only after it has been cleaned and is free from any materials which may constitute a health hazard. Use only cleansing agents which have been approved by Radleys to clean this device. To remove use:

Dyes	isopropyl alcohol
Construction materials	water containing tenside / isopropyl alcohol
Cosmetics	water containing tenside / isopropyl alcohol
Foodstuffs	water containing tenside
Fuels	water containing tenside

For materials which are not listed, please request information from Radleys. Wear the proper protective gloves during cleaning of the devices.

Electrical devices may not be placed in the cleansing agent for the purpose of cleaning.

Before using other than the recommended method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.

Associated standards and regulations

Construction in accordance with the foolwing safety standards

EN 61 010-1 UL 3101-1 CAN/CSA C22.2 (1010-1)

Construction in accordance with the foolwing EMC standards

EN 50 082-1 EN 55 014 EN EN60 555-2; -3 EN 50 081

Associated EU guidelines

EMC-guidelines: 89/336/EWG

Machine guidelines: 73/023/EWG

Technical data

design voltage:	VAC	230±10%
	or VAC	115±10%
design frequency:	Hz	50/60
motor:		EC-motor with ball bearings, without brushes
max. power input of motor	W	17
speed range:	l/min	50 bis 1100 infinitely variable setting silumin, GK AISi 12 polished
heating plate:		
diameter of heating plate:	mm	135
heating power:	W	600
max. surface temperature:	°C	300
max. surface temperature under error conditions	°C	370
temperature waviness of heating plate	K	±2,5
temperature setting:	°C	infinitely variable from ambient temperature
control accuracy of medium in ETS-D or contact- thermometer operation	K	±1 FUZZY-optimized
contact thermometer- connection permitted ambient temperature:	°C	DIN 12 878 Klasse 2 +5 bis +40
permitted relative humidity:	%	80
permitted duration of operation:	%	100
protection class acc. DIN 40 050:		IP 42
fuses:	F1/F2	10AT/250VId-Nr 2755400
	F3	1,025AT/250VId-Nr 2756300
contamination level:		2
overvoltage category:		II
protection class:		I (protective earth)
dimensions:	mm	160 x 280 x 95
weight:	kg	2,4

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List of spare parts Stirring Hotplate


Item	Designation	Item	Designation
2	bushing	46	plastic cover
3	radiating sheet	47	fixing plate
4	sealing washer	48	sealing washer
5	washer	102	ceramic stone
6	tapfite cheese-head screw	109	receptacle (gold-plated)
7	ceramic bushing	110	insulating hose
8	sleeve	2001	heating plate
9	foil keyboard	2002	EC-motor
10	housing-upper part	3001	poti
11	potentiometer shaft	3002	main plug
12	clamping bushing	3003	rocker switch
13	hood	3004	bushing plug
14	hexagonal nut	3005	rocker switch
15	cover for Pt 100	3006	relay
17	fixing clip	3007	fuse F1/F2
18	tapfite cheese-head screw	3008	fuse F3
19	toothed washer	5001	ball bearing
20	washer	5002	grooved ball bearing
21	O-ring (Perbunan)		
31	bottom part, coated		
32	plastite-screw		
33	BLP analog		
34	bushing		
35	magnet		
36	countersunk screw		
37	plastic cover		
38	tapfite countersunk screw		
39	cap		
40	screwed plug		
41	operating knob poti		
43	contact thermometer label		
44	contact plug		
45	tapfite cheese-head screw		

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Carousel Stirring Hotplates and Accessories

RR98072	Carousel Stirring Hotplate 230 volt & UK Plug
RR98072/EURO	Carousel Stirring Hotplate 230 volt & Euro Plug
RR98072/SWISS	Carousel Stirring Hotplate 240 volt & Swiss Plug
RR98072/USA	Carousel Stirring Hotplate 115 volt & US Plug
RR98072/JAP	Carousel Stirring Hotplate 110 volt & Japanese Plug
RR98073	Fuzzy Logic Digital Temperature Controller
RR98098	Spare Stainless Steel Sensor for Temperature Controller
RR98099	Glass Coated Sensor for Fuzzy Logic Temperature Controller
RR98075	Small Cross Shaped Stirring Bar, Pk 40
RR98091	Rare Earth Medium Cross Shaped Stirring Bar, Pk 20
RR98070	Large Octagonal Stirring Bar, Pk 20
RR98071	Small Octagonal Stirring Bar, Pk 40
RR98096	Rare Earth Small Elliptical Stirring Bar, Pk 40
RR98097	Rare Earth Large Elliptical Stirring Bar, Pk 20
RR98095	PTFE Stirring Bar Evaluation Kit, Box 18
RR98103	Stainless Steel Bath with Handle
RR98104	Stainless Steel Bath
RR98105	Special Sand Bath Filling
RR98094	Stirring Bar Retriever, 350mm
RR98106	Boss Head
RR98107	Holding Rod
RR98108	Support Rod
RR98109	Side Support Bracket
RR98110	Holding Fork and clamp
RR98111	Supporting Plate, 200mm diameter
RR98112	Thermomag Circulator
RR98100	Short Stainless Steel Sensor for Fuzzy Logic Temperature Controller
RR98102	1 metre Extension Cable for remote temperature sensing

**RADLEYS
DISCOVERY
TECHNOLOGIES**

Radleys Discovery Technologies Ltd
Shire Hill, Saffron Walden,
Essex, CB11 3AZ. United Kingdom.
Tel: +44-(0)1799-513320
Fax: +44-(0)1799-513283
Email: sales@radleys.co.uk
Web: www.radleys.com

IMPORTANT ELECTRICAL SAFETY NOTICE

Connection to the mains electrical supply

This equipment has been designed to one of the following safety classes. The appropriate warning should be observed.

Safety Class 1 THIS EQUIPMENT MUST BE EARTHED

Safety Class 2 THIS EQUIPMENT MUST NOT BE EARTHED

Before connecting the equipment to the mains electricity supply, please examine the information on the rating label affixed to the equipment.

Ensure that the supply is single phase alternating current (a.c) of the state frequency (Hz) with Neutral nominally at earth potential.

Check that the supply voltage is within the stated range.

The equipment rating label states the value of the input fuse fitted to the apparatus itself. Ensure that the plug or outlet circuit is fitted with an appropriate fuse of higher value.

For safety class 1

The wires in the mains lead are coloured in accordance to the following code.

Green and Yellow	-	Earth	(E)
Blue	-	Neutral	(N)
Brown	-	Live	(L)

Connect these wires to a non-reversible 3-pin plug as follows;

Green and Yellow wire to terminal marked: E (earth) or G (ground) or  or coloured Green or Green and Yellow

Blue wire to terminal marked: N (neutral) or common or coloured Blue

Brown wire to terminal marked: L (live) or Phase or coloured Brown

For safety class 2 - Double Insulated

The wires in the mains lead are coloured in accordance to the following code.

Blue	-	Neutral	(N)
Brown	-	Live	(L)

Connect these wires to a non-reversible 3-pin plug as appropriately directed above.

LIVE PARTS SHOULD NEVER BE EXPOSED UNLESS THE EQUIPMENT HAS BEEN SWITCHED OFF AND INSULATED FROM THE MAINS ELECTRICITY SUPPLY.